

# DermaSense

AI-Powered Skin Lesion Analysis  
& Bias Mitigation



Presented By:  
Romi Yosef  
Shir Molakandove  
Afik Haviv




# Phase 1: Automated Lesion Isolation

## Extracting Diagnostic Artifacts

The first stage of our pipeline focuses on isolating the biological lesion from its original light-skin context. This is critical to ensure that no "light-skin" features pollute the subsequent augmentation.

- **Library:** Leveraging the **RemBG** library for AI-powered background removal.
- **Masking:** Creating binary masks using *cv2.threshold* and *cv2.findNonZero* to locate precise coordinates.
- **Feature Preservation:** By isolating the mole, we preserve the original texture, asymmetry, and color distribution required for ResNet analysis.




# Phase 2: Synthetic Skin Synthesis

## Building the 'Digital Skin Bank'



To solve the data scarcity of dark skin tones, we utilize Stable Diffusion XL (SDXL) to generate a diverse library of synthetic backgrounds.

- **Medical Prompts:** Custom prompt engineering ensures the generation of *clinical-grade* textures, including pores, fine lines, and epidermal ridges.
- **Skin Fidelity:** We target the Fitzpatrick Scale Type VI to create the most underrepresented skin tone patches.
- **Precision:** Images are generated at 1024x1024 resolution and downsampled to match clinical input sizes.





# Phase 3: Intelligent Poisson Blending



## Seamless Image Fusion

The core of our augmentation is the Smart Blend function. It merges real malignant features with synthetic contexts without creating visual artifacts.

- Histogram Matching: We use `skimage.exposure` to synchronize the mole's lighting with the new dark skin background.
- Seamless Cloning: Implementing `cv2.seamlessClone` (Poisson Blending) to merge edges perfectly, maintaining diagnostic integrity.
- Result: A high-fidelity, medically valid dataset that represents cancer as it appears on dark skin.

# Phase 4 & 5: Biased vs. Diverse Training



## Establishing the Baseline

We first trained a ResNet50 model exclusively on original ISIC 2016 data. This confirmed our hypothesis: while accuracy was high on light skin, it achieved only 63.8% recall on dark skin, creating a dangerous clinical bias.



## Weighted Diverse Retraining

We retrained the model on our Hybrid Dataset. Crucially, we applied Class Weights in the loss function to force the model to prioritize malignant cases, regardless of the skin tone context.





## Model Refinement

Final optimization involved 25 epochs of fine-tuning with a low learning rate (0.0001), allowing the model to learn subtle melanoma markers on dark skin without forgetting light-skin features.

# Phase 6: Empirical Results & Impact

Biased Model DERMASENSE Model



## The Safety Win

By diversifying the dataset through generative synthesis, we achieved a **31.5% absolute increase** in malignant recall for minority populations. **DERMASENSE** successfully transforms AI from a biased tool into a life-saving assistant for everyone.

Thank you for your attention.

We are happy to take any questions.

